



SNTT1

Multipara Monitor

Key Features

- ✓ Waveforms on High-resolution Color TFT display.
- ✓ Graphical & tabular trends.
- ✓ In-built rechargeable battery.
- ✓ Interface facility with central monitoring system.
- ✓ Portable & easy to carry.
- ✓ Available in 7" Screen Size.



General

Dimensions	: 230mm X 200mm X 185mm (LXBXH)
Display	: 7" Color monitor
Resolution	: 800 x 480 pixels
User Interface	: Encoder

Power supply

Input voltage	: Adaptor 18V/1.2A
Battery type	: 12v/1.2AH (in-built rechargeable)
Charging time	: 80% charging in 12 hrs
Backup	: 2hrs on full charge

Environmental

Operating Temp	: -5 °C to +40 °C
Humidity	: up to 95 % non-condensing.

SpO2

Parameter	: SpO2%, PR, Bargraph, Perfusion index
Spo2 % Range	: 0 to 100%
Accuracy	: 70 to 100 +/-2%, 50 to 69+/-3%
Resolution	: 1 %
Indication	: Finger absent, Searching pulse
Pulse rate Range	: 10 – 250 BPM
Probe	: Adult, Pediatric, Neonate
Pitch variation change with oxygen saturation	

NIBP

Parameter	: Systolic, Diastolic and Mean
Method	: Oscillometric
Range	: 10 to 250
Time interval	: User selectable 3 min. to 250 min.
Modes	: Manual and Automatic
Cuff	: Single tube with Velcro.
Cuff deflation	: Automatic.
Patient Type	: Adult, Pediatric, Neonate

Settings

- o Limit setting for individual parameter
- o Facility to set default settings.

Keys:

Following hotkeys are provided for fast and easy operation

- o Alarm ON / OFF
- o QRS Beep ON / OFF
- o Take BP
- o Home
- o Power ON/OFF

knobs:

- o Encoder interface for browsing
- o Beep volume knob

*Due to our continuous product development program, the specification may change without prior notice.

NASAN MEDICAL ELECTRONICS PVT. LTD.

Plot No.147, "Rangali" Society, First Floor, Opp. Karve statue,

PUNE 411038, Maharashtra, INDIA Mobile: +918237032389

Email: sales@nasanmedical.com Website: www.nasanmedical.com